

Open Wide® Mouth Prop for those with Special Needs in Oral Hygiene



**Patient-friendly Mouth Prop has
a long handle for better control**



**The OWD is the perfect adjunct to
the Surround® Toothbrush**

- *Wedge-shape accommodates most patients*
- *Foam handle provides comfort and control*
- *Rigid wooden core adds stability and strength*
- *Handle-end can also be used for cheek retraction*



The **Open Wide® Mouth Prop** (aka mouth rest) from Specialized Care Co*. USA was designed for those who are unable to cooperate or have difficulty with motor control in keeping their mouths open during oral care. Perfect companion product to the **Surround® Toothbrush**, the OWD makes oral hygiene safer and less stressful for both parent /carer and patient.

Simple to use and effective in maintaining oral access, the OWD has a graduated wedge design with padded handle for stability & comfort; long handle is easy to grip and maintain control during brushing. Made from a special high-density foam that is durable and difficult to bite through, the OWD has a rigid wooden core for added stability, yet is comfortable and gentle in the mouth.



**40%
More
Foam!**

The OWD can be cleaned in hot water or dishwasher, reused repeatedly for the same person and replaced as needed. (*Assess integrity of the outer foam after each use, "dimpling" from the teeth may occur as the foam depresses when bitten*)

For those with more aggressive bites a "thicker" version is available with 40% more foam for added strength. (*OWD-T*)

Long favoured by Special Needs Dentists, the OWD is the preferred mouth prop to facilitate oral care in Nursing Homes, hospitals and home health settings.