

Stay N Place® Airway Positioners


Our Airway Positioners Help When "Chin Up" is a Challenge.



Maintaining proper positioning for procedures can be challenging. Our **Stay N Place® Airway Positioners** support the upper torso, making "chin up" easier. The **Stay N Place® Airway Positioner (code: SP-AP VIN)** above, has two head and neck cradles and a wedge-shaped pillow to lift the chest area. These vinyl-covered pillows can be repositioned to accommodate a range of patients. It can be used with or without the head strap. Available in dark teal or cream vinyl.



The foam **Stay N Place® Airway Positioner (code: SP-AP FM)** provides consistent support (less compressibility), which is important for sedation as well as routine operative dentistry. The surface is smooth to the touch, hygienic and is easily cleaned.



Stay N Place® Airway Positioners

Stay N Place® Airway Positioner-Vinyl

SP-AP VINdt (dark teal)

SP-AP VINc (cream)

Stay N Place® Airway Positioner - Foam

SP-AP FM (soothing blue)

100 Protective Sheaths - SH-100/AP

Creating dental products that help you provide safe, comfortable care.

www.specializedcare.com

International Distributor: Mobile Medical Systems International /02 9252 4410/info@mobilemedical.com.au

Table of Contents

February 2018

Positioning Products

Stay N Place® Airway Positioners.....	Page 1
Stay N Place® Head and Neck Rest.....	Page 2
Stay N Place® Lap Cushion/Stay N Place® Lap Board.....	Page 3
Stay N Place® Booster Seat/Rainbow® Toe Guard.....	Page 4
Stay N Place® Chair Cushions.....	Page 5
Stay N Place® Pillows.....	Pages 6 & 7

Home Care Products

Surround® Toothbrushes.....	Pages 8 & 9
Open Wide® Reusable Mouth Prop.....	Page 10
Open Wide® Mouth Rest	Page 11
Oral Care Kits & Home Training Videos.....	Page 12

Stabilizing Products

Professional/Parent Education, D-Termined Program.....	Page 13
Robot Arm™ Stabilizers.....	Page 14
Rainbow® Boards.....	Page 16
Rainbow® Wraps.....	Page 17
Rainbow® Accessories.....	Page 18
Patient Education(Stabilizing Products).....	Page 19

We accept: Visa  Mastercard  American Express  Discover 



Specialized Care Co, Inc.
....from Design to Delivery

Specialized Care Company's mission is to design, manufacture and market dental products that serve the unique needs of dentists, hygienists, health care providers, caregivers and parents who provide dental care for young children and people with special health care needs.

We understand that we have an important responsibility to our customers. Our products are tested in a realistic clinical environment. We listen to our customers so that we may continually improve existing products and provide needed new ones.

Orders must be prepaid with Money Order, Credit Card or check. TERMS are Net 30 days for approved facilities. Federal ID # 02-051-2077

Stay N Place® Head & Neck Rest

Novel Backboard Supports the Patient in Her Own Chair.

When it's best for the patient to remain in her own chair during dental treatment, the **Stay N Place® Head & Neck Rest** provides support as her head is tilted back. Result? Dental care can be safely delivered without transferring.

The Stay N Place Head & Neck Rest is positioned behind the patient and held in place with straps secured to the chair. The bottom of the device can be brought forward (away from the back of the chair) to achieve more of an angle.

This low-tech device is a real multi-tasker. It can be used for in-office as well as off-site care. All of the pillows can be repositioned for patient comfort, and any of the pillows can be removed as needed. Clean with disinfecting agent and wipe dry with a soft cloth.



Stay N Place® Head & Neck Rest	Dark Teal	Cream
Stay N Place® Head & Neck Rest	SPHNRdt	SPHNRc

Creating dental products that help you provide safe, comfortable care.

www.specializedcare.com

International Distributor: Mobile Medical Systems International /02 9252 4410/info@mobilemedical.com.au


Stay N Place® Lap Cushion

The Lap Cushion Supports Young Patients During An Infant Oral Exam.

Many professionals use the knee-to-knee technique for infant oral exams. The **Stay N Place® Lap Cushion** makes the exam more secure for the baby, more comfortable for the adults and more hygienic for all. The cushion is held under the patient as he is tilted into the lap of the professional, thus increasing the baby's sense of protection. The dentist has proper access for the exam and the caregiver is assured. Parents appreciate the cushy feel and the secure side rails. Clean with disinfecting agent and wipe dry with a soft cloth.



The Stay N Place Lap Cushion can be held behind the child as he is tilted back. Mom is in sight, which is reassuring to the baby.



Flexible, comfy and durable, our latest infant lap cushion has it all!

Stay N Place® Lap Cushion

Stay N Place® Lap Cushion	SP-FLEX	24"L x 17"W
100 Protective Sheaths - SH-100/LAP		12" x 8" x 30" gusseted


Stay N Place® Lap Board

Our Lap Board Provides Another Choice for Knee-to-Knee Exams.

Our **Stay N Place® Lap Board** also assists with infant oral exams by providing a hygienic surface to support the baby. The difference is the construction. The Stay N Place Lap Board is made with a rigid platform, layered with foam and upholstered in soft, durable vinyl. Raised sides help keep baby safely on the Lap Board.



The Stay N Place Lap Board provides a safe, rigid surface for the child while allowing distance between the adults.



Stay N Place® Lap Board

	Dark Teal	Cream
Stay N Place® Lap Board (24"L x 16"W)	SP-LAP BDdt	SP-LAP BDC
100 Protective Sheaths - SH-100/LAP	12" x 8" x 30" gusseted	

Creating dental products that help you provide safe, comfortable care.

www.specializedcare.com

International Distributor: Mobile Medical Systems International /02 9252 4410/info@mobilemedical.com.au

Stay N Place® Booster Seat

The Booster Seat Saves Your Back and Provides Comfort for Young Patients.



Removable lumbar support

Do you find yourself “scooting” patients to the top of the chair for treatment, only to have them slide back down when the chair is brought upright? Daily repetition of that scenario wreaks havoc on your back, causes lost time for the practice and is unsettling for the patient.

The Stay N Place® Booster Seat positions your patient 5” higher, providing better access and increasing efficiency. This simple addition to the dental chair is more ergonomically correct for you and more comfortable for your patient. It now comes with a removable lumbar support pillow.

The Stay N Place Booster Seat can be cleaned with the same disinfectant used for vinyl-covered dental chairs.

The Stay N Place Booster Seat can be used with or without the lumbar support.

Stay N Place® Booster Seat	Dark Teal	Cream
Stay N Place® Booster Seat	SPBSdt	SPBSc

Important: It is the clinician's responsibility to assure the appropriateness of any medical device for a given patient. Be sure to understand any applicable regulations.

Rainbow® Toe Guard

A Cleverly Designed Pillow Keeps Patient from Interfering with Monitors.



Left: The patient can dislodge the monitoring device.



Right: The Rainbow Toe Guard is an effective yet comfortable way to solve the problem.

The Rainbow® Toe Guard prevents the patient from using one foot to kick a monitoring device off the other foot. The specially designed pillow is placed between the feet and secured with ankle straps. The belt is used around the pillow and stabilizing board or around the dental chair. This is a reliable way to keep the patient monitored and comfortable.

Rainbow® Toe Guard	Dark Teal	Cream
Rainbow® Toe Guard - One size	RTGdt	RTGc

Creating dental products that help you provide safe, comfortable care. www.specializedcare.com

International Distributor: Mobile Medical Systems International /02 9252 4410/info@mobilemedical.com.au

Stay N Place® Chair Cushion

This Cushion Provides Physical and Emotional Comfort.



Both this chair cushion and the toddler version come with Roll-N-Store Straps for easy carrying and storage.

When body shape and chair shape are not a good match, or when patients become fearful at the sight of a dental chair, the **Stay N Place® Chair Cushion** can help. It presents a welcoming image to help calm patient anxiety. Crushed foam displaces as the patient settles in, filling gaps in the chair and providing physical comfort and a secure feeling. The Stay N Place Chair Cushion curves in at the neck, allowing your patient's head to tilt back and the chest to be lifted, while leaving you room to maneuver.

A non-slip bottom helps the cushion stay where you place it on the dental chair. Available in soft, yet durable vinyl with a choice of dark teal or cream.

Stay N Place® Chair Cushion	Dark Teal	Cream
Stay N Place® Chair Cushion	SPCCdt	SPCCc

Important: It is the clinician's responsibility to assure the appropriateness of any medical device for a given patient. Be sure to understand any applicable regulations.

Stay N Place® Chair Cushion For Toddlers

The Infant Chair Cushion Makes it Easier to Accommodate Toddlers.

An over-sized dental chair can feel insecure for toddlers, and it can cause you to bend and twist in order to provide treatment. The **Stay N Place® Chair Cushion for Toddlers** creates a safe, cozy lift for the patient, presents a professional image for parents, and helps you maintain better posture.

The chair cushion is like a bean bag, with soft foam inside. Air escapes through side vents as the child is placed upon it. This settles the cushion into the curves of the dental chair while securely nestling the youngster. It has a non-slip bottom, and comes with a safety belt for added security.

Designed for dentistry, the Stay N Place Chair Cushion has a soft, durable vinyl cover. Clean with disinfecting agent and wipe dry with a soft cloth.



Before: The dentist is bending and twisting to do her job



After: The patient is at a more ergonomic height for dental treatment.

Stay N Place® Chair Cushion/ Toddlers	Dark Teal	Cream
Stay N Place® Chair Cushion/Toddlers	SPCC INFdt	SPCC INFc

The Stay N Place Chair Cushion for Toddlers brings the patient to a more ergonomically correct height.

Creating dental products that help you provide safe, comfortable care.

www.specializedcare.com

International Distributor: Mobile Medical Systems International /02 9252 4410/info@mobilemedical.com.au

Stay N Place® Pillows

Specially Designed Pillows Support the “Chin Up” Position.

No chair can accommodate all patients equally well. If the patient falls short of the chair’s headrest, he or she may not be properly supported. **Stay N Place® Pillows** provide the support you need to do your job effectively and the comfort the patient needs to sit without fidgeting.



The Neck Pillow supports the neck, helps keep the head tilted back, and can be used as a shoulder roll.

Pediatric - 3.5” H x 3.5” W x 10.5”L
 Adult - 4.5” H x 4.5” W x 12.5”L



The sides of our Head Cradle push up into a cradling position when the patient puts pressure on the center section.

Pediatric - 3” H x 3.5” W x 9.5”L
 Adult - 4” H x 4” W x 12”L



The Knee Pillow helps to center the patient and relieves pressure on the lower back.

Pediatric - 4.5” H x 4.5” W x 12.5”L
 Adult - 6” H x 6” W x 20”L



All pillows are available in dark teal or cream vinyl.

A. Stay N Place® Neck Pillow	Dark Teal	Cream
Neck Pillow - Pediatric	SPNP-Pdt	SPNP-Pc
Neck Pillow - Adult	SPNP-Adt	SPNP-Ac
B. Stay N Place® Head Cradle	Dark Teal	Cream
Head Cradle - Pediatric	SPHC-Pdt	SPHC-Pc
Head Cradle - Adult	SPHC-Adt	SPHC-Ac
C. Stay N Place® Knee Pillow	Dark Teal	Cream
Knee Pillow - Pediatric	SPKP-Pdt	SPKP-Pc
Knee Pillow - Adult	SPKP-Adt	SPKP-Ac

Creating dental products that help you provide safe, comfortable care.

www.specializedcare.com

International Distributor: Mobile Medical Systems International /02 9252 4410/info@mobilemedical.com.au

Stay N Place® Pillows

Support in the Right Places Facilitates Good Clinical Outcomes.

All our pillows have a strap that attaches behind the dental chair to keep the pillow in place, even when the patient sits forward. The vinyl pillow is soft to the touch, resists fluids and is easily cleaned using a surface disinfectant. Pillows are available in dark teal or cream.



D Head Pillow

The Head Pillow is designed for comfort and gives your patient a sense of security.

One size - 4" H x 11" W x 14" L



E Conforming Pillow

The Conforming Pillow takes the shape of the patient's head and neck and has a center cutaway to accommodate pulled-back hairstyles.

One size - 4" H x 7" W x 9" L

Tip of the day:

When using our website,
sign into your account
BEFORE shopping to see
codes and pricing
reserved for professionals.



F Plastic Sheaths

Plastic sheaths have a hole to allow the pillow strap to pass through and will fit all pillow configurations.

D. Stay N Place® Head Pillow	Dark Teal	Cream
Head Pillow	SPHdt	SPHc
E. Stay N Place® Conforming Pillow	Dark Teal	Cream
Conforming Pillow	SPCPdt	SPCPc
F. Stay N Place® Plastic Sheaths		
Package of 100 Sheaths for Pillows	SH-100	

Replacement vinyl covers are available for all pillows.
Please inquire.

Creating dental products that help you provide safe, comfortable care.

www.specializedcare.com

International Distributor: Mobile Medical Systems International /02 9252 4410/info@mobilemedical.com.au

Surround[®] Toothbrush

This Unique Design Cleans All Surfaces of the Teeth at the Same Time.



The Surround[®] Toothbrush is a great choice for patients with low tolerance for daily toothbrushing. It has three bristle pads to clean all surfaces of the teeth at the same time. Whether the user is brushing her own teeth or someone else is doing the brushing, the Surround Toothbrush can make the process easier because more tooth surface can be cleaned with each stroke.



A soft, rubber-like material covers the length of the toothbrush to help protect the oral environment if the person bites down unexpectedly, or falls while using the brush. The soft, end-rounded bristles and heart-shaped, compact head make brushing safe and comfortable.



Also available in a toddler size (up to age 5). See next page.

Latex-Free *

All toothbrushes come with a reusable plastic cover.

Regular	Code
Single Toothbrush	STB-1
3 Pack	STB-3
6 Pack	STB-6
12 Pack - Professional Sale Only	
72 Pack - Professional Sale Only	
144 Pack - Professional Sale Only	

* Not made with natural rubber latex

Tip of the day:
 When using our website, sign into your account **BEFORE** shopping to see codes and pricing reserved for professionals.

Creating dental products that help you provide safe, comfortable care. www.specializedcare.com

International Distributor: Mobile Medical Systems International /02 9252 4410/info@mobilemedical.com.au

Surround® Toddler Toothbrush

This Toothbrush Helps Teach Young Patients How to Brush.



The **Surround® Toddler Toothbrush** is designed with 3-bristle pads arranged at the correct angles for proper cleaning of teeth and gums. As the child attempts independent brushing, more surface is cleaned with each brush stroke because the bristles are in contact with lingual, facial and occlusal surfaces at the same time.

While made for toddlers, this brush is great for people of all ages who wear braces. Shorter bristles get in and around the appliance.

Our Home Care Products Are Also Available On Amazon.Com

Surround® Toothbrush • Surround® Toddler Toothbrush

Open Wide® Mouth Rest - Regular or Thick

Kit with "How to Brush the Teeth of Another Person " DVD

Kit sold on Amazon contains 2 Open Wide Mouth Rests, 1 Surround Toothbrush, and an Instructional DVD

Toddler	Code
Single Toothbrush	STBC-1
3 Pack	STBC-3
6 Pack -	STBC-6
12 Pack - Professional Sale Only	
72 Pack - Professional Sale Only	
144 Pack - Professional Sale Only	

* Not made with natural rubber latex

Tip of the day:

When using our website,
sign into your account
BEFORE shopping to see
codes and pricing
reserved for professionals.

Creating dental products that help you provide safe, comfortable care.

www.specializedcare.com

International Distributor: Mobile Medical Systems International /02 9252 4410/info@mobilemedical.com.au

Open Wide® Reusable Mouth Prop

These Mouth Props Give the Clinician Control for Proper Care.



The *Wraparound* handle allows the clinician to control the prop and the patient's mandible at the same time.

Do you use a mouth prop for patients who have difficulty keeping the mouth open reliably? Those same patients are the ones most likely to attempt to dislodge the prop. Our unique **Wraparound** mouth prop allows the clinician to maintain control with a handle that wraps to the outside of the cheek so it will not block oral access. The mouth prop is designed with a perfectly placed hole to keep a saliva ejector where it belongs. The prop is made of plastic that is overmolded with a crush-resistant latex-free* material.



The **Tab** mouth prop is streamlined for procedures that need greater clearance. It has a hole in the tab to allow floss to be threaded for easy retrieval.



Both styles can be autoclaved repeatedly without getting "gummy."



Shown are sizes Large (green), Medium (blue) and Small (purple) with **Wraparound** handle.

Open Wide® Mouth Prop with Wraparound Handle

OWR-2S	2 Small Mouth Props/WRAPAROUND
OWR-2M	2 Medium Mouth Props/WRAPAROUND
OWR-2L	2 Large Mouth Props/WRAPAROUND
OWR-S/M	1 each Small & Medium Mouth Prop/WRAPAROUND
OWR-M/L	1 each Medium & Large Mouth Prop/WRAPAROUND
OWR/3 SZS	1 each Small, Medium & Large Mouth Prop/WRAPAROUND



Shown are sizes Large (green), Medium (blue) and Small (purple) with **Tab** handle.

Open Wide® Mouth Prop with Tab Handle

OWT-2S	2 Small Mouth Props/TAB
OWT-2M	2 Medium Mouth Props/TAB
OWT-2L	2 Large Mouth Props/TAB
OWT-S/M	1 each Small & Medium Mouth Prop/TAB
OWT-M/L	1 each Medium & Large Mouth Prop/TAB
OWT/3 SZS	1 each Small, Medium & Large Mouth Prop/TAB

* Not made with natural rubber latex

Important: It is the clinician's responsibility to assure the appropriateness of any medical device for a given patient. Be sure to understand any applicable regulations.

Creating dental products that help you provide safe, comfortable care.

www.specializedcare.com

International Distributor: Mobile Medical Systems International /02 9252 4410/info@mobilemedical.com.au

Open Wide® Mouth Rest

A Reliably Open Mouth is Important for In-Office and At-Home Care.

The **Open Wide® Mouth Rest** is a helpful tool for those who lack the cooperation or motor control to keep the mouth reliably open. It is made from special high-density foam with a rigid core to add stability. It is sturdy enough to resist biting pressure, but has enough “give” to be comfortable for the person receiving care. The mouth rest has a long handle for safe, easy control. It is considered a disposable item in the office but can also be cleaned and given to a trained caregiver for daily oral hygiene.

In-Office



At-Home



Important: Consult the appropriate authority in your state or facility for regulations that may pertain to the use of mouth rests.

Open Wide® Mouth Rest

Description

OWD-4	4 pack (Regular)
OWD-24 (in plastic storage container)	24 pack (Regular)
OWD-50	50 pack (Regular)
OWD-100	100 pack (Regular)

Open Wide® Mouth Rest Thick

Description

OWDT-4	4 pack (Thick)
OWDT-24 (in plastic storage container)	24 pack (Thick)
OWDT-50	50 pack (Thick)
OWDT-100	100 pack (Thick)

U.S. Patent No. 5,735,691



The Thick Mouth Rest has 40% more foam to handle aggressive biting pressure.

Creating dental products that help you provide safe, comfortable care.

www.specializedcare.com

International Distributor: Mobile Medical Systems International /02 9252 4410/info@mobilemedical.com.au

Robot Arm™ Stabilizers

Help Teach Young Patients How To Sit Safely for Dentistry



These arm and leg stabilizers control kicking and grabbing movements that may pose a danger to both patient and practitioner.

The **Robot Arm™ Stabilizers** are a perfect tool to teach young patients how to keep hands out of harm's way. They start out as flat sheets and are formed into cylinders around arms or legs to prevent grabbing, kicking and bending movements. The Velcro® material makes the arm stabilizers easy to apply and remove during the appointment as a patient learns. They are colorful and fun, appealing to parents and patients alike.

The Robot Arm Stabilizers are just one in a series of products that help you provide safe, comfortable care to children and people with special needs. Please see the following pages for information about other products that help prevent sudden, unsafe patient movement.

Important: It is the clinician's responsibility to assure the appropriateness of any medical device for a given patient. Be sure to understand all applicable regulations regarding the use of protective stabilization. Never leave patient unattended.

Robot Arm™ Stabilizers (packaged in pairs)

Code	Size (Circumference x Height)	Usage Guide
RAS-SM	10" x 8" (up to 9" arms/legs)	Ages 1-4
RAS-MED	12" x 10" (up to 11" arms/legs)	Ages 5-8
RAS-LG	14" x 12" (up to 13" arms/legs)	Ages 9-12
RAS-XL	16" x 18" (up to 15" arms/legs)	Teens
RAS-XXL	20" x 18" (up to 19" arms/legs)	Adults

Velcro® is a registered trademark of Velcro BVBA.

Creating dental products that help you provide safe, comfortable care.

www.specializedcare.com

International Distributor: Mobile Medical Systems International /02 9252 4410/info@mobile-

www.specializedcare.com

Rainbow® Stabilizing System

Protective Stabilization Helps Keep Everyone Safer.

Whether you are treating a very young child, or an older person with behavioral limitations, an unavoidable reality is that some patients may become unexpectedly active during a procedure. Since dentistry involves the use of sharp instruments in a small, highly vascularized area, you owe it to your patients and to your practice to use appropriate safety equipment.

The Rainbow® Stabilizing System is a state-of-the-art protective device consisting of a colorful, breathable mesh wrap and vinyl-covered, padded board to increase acceptance by both patient and caregiver.



Our young model helps to demonstrate the use of the small Rainbow Stabilizing System. Reagan is 3 years old, 3 feet and about 31 pounds.

See pages 16-18 for more details on our Rainbow® Stabilizing System and accessories.

Creating dental products that help you provide safe, comfortable care.

www.specializedcare.com

International Distributor: Mobile Medical Systems International /02 9252 4410/info@mobilemedical.com.au

Rainbow® Board

The Rainbow® Board Helps You Put Safety First During Dental Procedures.

The **Rainbow® Board** provides the foundation for the stabilizing system. It is available as a **Hinged Board** or a **Straight Board**. The Hinged Board bends with the dental chair and allows the patient to be moved to an upright position quickly, which is helpful should the patient have a medical emergency during a procedure. We offer the Straight Board for practitioners whose equipment is better suited to a rigid board, such as a bench-style chair. Rainbow Boards are used in conjunction with Rainbow Wraps (next page).

The vinyl covered **Rainbow Board** is padded for patient comfort.



The helpful **Rainbow Head Tilt** keeps the upper torso lifted and the head in the “chin up” position.

The multi-use **Rainbow Knee Support** relieves pressure on the lower back. Both are sold as accessories to the board.

Rainbow® Head Tilt & Knee Support

- Rainbow® Head Tilt - Small/Medium
- Rainbow® Head Tilt - Large
- Rainbow® Knee Support - One Size

Dark Teal

- RHT-S/Mdt
- RHT-Ldt
- RKSdt

Cream

- RHT-S/Mc
- RHT-Lc
- RKSc



The Hinged Boards fold for easy storage.



The Large Hinged Rainbow Board articulates in the center and at the bottom to conform to the contours of the dental chair.



Rainbow® Boards

Height	Width	Rainbow® Board	Dark Teal	Cream
2 Feet	10"	Infant Board	RBldt	RBlc
3 Feet	10"	Small Board	RBSdt	RBSc
4 Feet	11"	Medium Hinged Board	RBHMdt	RBHMdc
4 Feet	11"	Medium Straight Board	RBMDt	RBMc
5 Feet	15"	Large Hinged Board	RBHLdt	RBHLc
5 Feet	15"	Large Straight Board	RBLdt	RBLc

Clean with a disinfecting agent and wipe dry with a soft cloth.

Creating dental products that help you provide safe, comfortable care. www.specializedcare.com

International Distributor: Mobile Medical Systems International /02 9252 4410/info@mobilemedical.com.au

Rainbow® Wrap

Our Unique Wrap is Designed Specifically for Dentistry.

The **Rainbow® Wrap** keeps arms and legs out of harm's way while sharp instruments are in use. Everything about it has been designed for the needs of a dental office. The breathable mesh avoids overheating that often accompanies patient agitation. An opening in the wrap allows you to conveniently use monitoring devices. The wrap has a friendly look, so it is better accepted by patients and parents. But it doesn't stop there. The wrap is made of a moisture-resistant, washable fabric - so even laundering is made simple.



The knee belt can be attached to the wrap with snaps for more efficient preparation.

New! Tapered shoulders for a more secure fit.

Important: It is the clinician's responsibility to assure the appropriateness of any medical device for a given patient. Be sure to understand all applicable regulations regarding the use of protective stabilization. Use appropriate monitoring, and never leave patient unattended.

What Size to Buy?

Up to this approx. weight/age:

Recommendation

26 pounds.....	(Infants).....	Infant Board - Infant Wrap
35 pounds.....	(1-4 years old).....	Small Board - Small Wrap
60 pounds.....	(5-8 years old).....	Medium Board - Medium Wrap
145 pounds.....	(9-12 years old).....	Large Board - Large Wrap
185 pounds.....	(Teen/Adults).....	Large Board - X-Large Wrap
250 pounds.....	(Teen/Adults).....	Large Board - XX-Large Wrap

Recommended Board and Wrap sizes, based on approximate age and weight of patient.



New! Our latest Wristlets are cushioned for patient comfort.

Rainbow® Wristlets

RWC10/PR	One pair (10")
RWC13/PR	One pair (13")

Rainbow® Wraps with easy-clean back panel

Code	Rainbow® Wrap	Length/Width	Wrap Comes With.....
RWlp	Infant Wrap	20"/43"	Wristlets, Laundry Panel, Storage Bag*
RWSp	Small Wrap	25"/35"	Wristlets, Laundry Panel, Storage Bag*
RWMP	Medium Wrap	32"/44"	Knee Belt, Wristlets, Laundry Panel, Storage Bag
RWLp	Large Wrap	34"/49"	Knee Belt, Wristlets, Laundry Panel, Storage Bag
RWXLp	X-Large Wrap	44"/55"	Knee Belt, Wristlets, Laundry Panel, Storage Bag
RWXXLp	XX-Large Wrap	44"/70"	Knee Belt, Wristlets, Laundry Panel, Storage Bag

* Knee Belt Not Required for Small or Infant Wraps

Creating dental products that help you provide safe, comfortable care.

www.specializedcare.com

International Distributor: Mobile Medical Systems International /02 9252 4410/info@mobilemedical.com.au

Accessories for Your Rainbow® System

Accessories Help Customize the Rainbow® System for Your Practice.



The Rainbow® Head Stabilizer is a comprehensive head control device that minimizes both up-and-down and side-to-side movement. It attaches to the Rainbow Board with screws for reliable stability. Adjustable side panels tilt in at a slight angle and accommodate a wide range of patients. The thick cushions on the side panels allow for a snug fit while maintaining patient comfort. The Head Stabilizer must be used in conjunction with the Rainbow Board.

Rainbow® Head Stabilizer	Dark Teal	Cream
Rainbow® Head Stabilizer for Infant Board	RHlIdt	RHlIc
Rainbow® Head Stabilizer for S/M Boards	RHlCdT	RHlCc
Rainbow® Head Stabilizer for Large Boards	RHlAdt	RHlAc
Replacement Head Strap - RHS-2		

The Rainbow® Triangle Pillow fills the gap in the dental chair where the seat meets the back. It is covered with a non-slip material, so it stays where you put it. Our customers report it also makes a handy knee pillow.



Rainbow Triangle Pillow close-up (above) and in use (right).

Rainbow® Triangle Pillow	Dark Teal	Cream
Rainbow Triangle Pillow	RTPdt	RTPc



Rainbow Knee Belt



Rainbow Wide Safety Belt

Both the **Rainbow® Knee Belt** and the **Rainbow® Wide Safety Belt** control leg movements that can be a danger to patient and staff. The latter is longer to be wrapped around a dental chair.

Rainbow® Knee Belt and Rainbow® Wide Safety Belt		
RKB-Mp	Medium Knee Belt	(6"W by 40"L)
RKB-Lp	Large Knee Belt	(9"W by 43"L)
RKB-XLp	X-Large Knee Belt	(9"W by 58"L)
RWSB	Rainbow® Wide Safety Belt	(4"W by 60"L)

Rainbow® Storage Bags - Canvas-like fabric

- RST-S Rainbow® Storage Bag for Small or Infant System
- RST-M Rainbow® Storage Bag for Medium System
- RST-L Rainbow® Storage Bag for Large System
- RST-H Rainbow® Storage Bag for Hinged Systems



Rainbow® Storage Bags store your board, wrap, belts, etc. They are made of a tough canvas-like fabric to extend the life of your investment, and come with a carry strap.

Creating dental products that help you provide safe, comfortable care.

www.specializedcare.com

International Distributor: Mobile Medical Systems International /02 9252 4410/info@mobilemedical.com.au

Patient Education

These Educational Tools Help You Explain Stabilization.

The **Protective Stabilization Model** is a miniature stabilizing system complete with “patient”. It helps you demonstrate to caregivers how a patient is safely swaddled during treatment. It will become a key part of your informed consent procedure especially when there is a language barrier. The model measures 11” tall by 6.5” wide, and can also be used to help patients know what to expect.



Protective Stabilization Model

PSM	Protective Stabilization Model
HB-Pkg	Educational Package 1 PSM and 3 <i>Riley and the Hugging Blanket</i> Books
HBS-Pkg	Educational Package (Spanish) 1 PSM and 3 <i>Raffi Y La Cobija Abrazadora</i> Books



Riley and the Hugging Blanket is a beautifully illustrated book about a patient who gets upset at the dental office, despite trying to be brave. The “Hugging Blanket” is explained to Riley, who accepts its use. Riley learns to relax, and soon the Blanket is no longer needed. Riley is very proud!

The teaching book is available in English and in Spanish. It can be used by parents at home to acclimate the child to the concept of stabilization prior to an appointment.

Riley and the Hugging Blanket Children's Book

HBK-1	<i>Riley and the Hugging Blanket</i> Book - Single Copy
HBK-3	<i>Riley and the Hugging Blanket</i> Book - 3 Pack

Raffi Y La Cobija Abrazadora Children's Book

HBK-S1	<i>Raffi Y La Cobija Abrazadora</i> - Single Copy
HBK-S3	<i>Raffi Y La Cobija Abrazadora</i> - 3 Pack



Please Note: We make some of our products available directly to at-home users on our website. As a professional, we ask you to sign in to your account *before shopping* so that you will see items and pricing reserved for professionals.

Creating dental products that help you provide safe, comfortable care.

www.specializedcare.com

International Distributor: Mobile Medical Systems International /02 9252 4410/info@mobilemedical.com.au

Return Service Requested



We accept: Visa  Mastercard  American Express  Discover 

Orders must be prepaid with Money Order, Credit Card or check. TERMS are Net 30 days for approved facilities. Federal ID # 02-051-2077

Solutions for Your Everyday Challenges

Patient won't keep the mouth open	<ul style="list-style-type: none"> • Mouth Props..... Page 10 • Mouth Rests.....Page 11
Patient does not maintain "chin up" position	<ul style="list-style-type: none"> • Airway Positioners.....Page 1 • Dental Chair Pillows..... Pages 6 & 7 • Head Tilt and Knee Support.....Page 16
Patient lacks cooperation skills	<ul style="list-style-type: none"> • Stabilizing System..... Pages 16 & 17 • Robot Arm™ Stabilizers..... Page 14
Patient does not fit well in the dental chair	<ul style="list-style-type: none"> • Booster Seat..... Page 4 • Chair Cushions..... Page 5 • Dental Chair Pillows..... Pages 6 & 7
Caregiver or patient needs help understanding stabilizers	<ul style="list-style-type: none"> • Protective Stabilization DVD.....Page 13 • Protective Stabilization Model..... Page 19 • "Riley and the Hugging Blanket" Book..... Page 19
Infant patient requires Knee-to-Knee Exam	<ul style="list-style-type: none"> • Lap Cushion.....Page 3 • Lap Board.....Page 3
Patient tries to kick pulse oximeter off toe	<ul style="list-style-type: none"> • Toe Guard.....Page 4
Not safe to transfer the patient from wheelchair	<ul style="list-style-type: none"> • Head & Neck Rest..... Page 2
Patient does not cooperate for daily oral hygiene	<ul style="list-style-type: none"> • Mouth Rests.....Page 11 • Toothbrushes.....Pages 8 & 9

Creating dental products that help you provide safe, comfortable care. www.specializedcare.com

International Distributor: Mobile Medical Systems International /02 9252 4410/info@mobilemedical.com.au